ARAHAN UNTUK KEGUNAAN

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FFP3V R D SAFETY RESPIRATOR

UK Technology. Manufactured In Malaysia !

UK Protection Factor
APF20 x WEL

Conforms to EN149

**SIRIM
TESTED
CERTIFICATE
ISSUANCE
IN PROGRESS ***

FARZIM MEDICARE Sdn Bhd
T : +603-2785 6842 F : +603-2785 6942
farzinhq@farzimmcare.com.my
Website : www.farzimmcare.com.my
No 308A, Level 1, Lorong Selangor, Taman Melawati,
53100 KUALA LUMPUR

Global brands of quality medical devices.
Previously known as
First Aid Red Zone Instrument Malaysia (F-A-R-Z-I-M)
Medical Device Authority (MDA) Reg. No. : KP 570 828 7315
ROC Reg. No. : 200901022838 (865937-P)

**FARZIM
MEDICARE**

...its for the best, ...only from the best !*

*email inquiries to SIRIM QAS International Sdn Bhd at aishahs@sirim.my

PRODUCT INSERT INSTRUCTION FOR USE

INSTRUCTIONS FOR USE

This NAFAShield 99 FFP3 V R D Safety Respirator is certified in accordance with EN 149. Please refer to our Product Catalogue for more information on the product certification and performance.

GUIDELINES FOR THE SELECTION AND USE MAY VARY ACCORDING TO INDIVIDUAL COUNTRIES. ALWAYS CHECK PRIOR TO USE TO ENSURE THAT YOUR SELECTION FOR THE INTENDED APPLICATIONS MEETS NATIONAL REQUIREMENTS.

APPLICATION TABLE

NAFAShield 99 FFP3 V R D Safety Respirator is suitable for the following applications:

Filter Type : EN149-FFP3 Colour : WHITE	Application Type : Particulate Filter
Filter Class : FFP3 Minimum 99% Efficiency @0.3 micron	Main Application : For use against hazardous particulates both solid and liquid particles. Dust, mist, fumes of all types and micro-organisms e.g. bacteria and viruses.

CHECK PRIOR TO USE

- Prior to use ensure that the mask has been stored in accordance with the recommended storage instructions below and are still within their stated expiry date, the detail is displayed on the carton label.
- Remove all packaging and inspect both the filter and the mask for damage. If in doubt replace any suspected defective item.
- Ensure that the intended application and atmosphere has been assessed by a competent person and meets the application limitations for the NAFAShield 99 FFP3 V R D Safety Respirator (see application table and the additional warnings & limitations).

FITTING INSTRUCTIONS

- Ensure the checks prior to use above have been conducted before the use of each mask.
- Remove eyewear. If facial hair protrudes beyond the border of the mask an effective seal will not be achieved and a powered respirator will be required.
- Grasp the respirator with one hand positioning your thumb on the inside of the mask. (Fig 1a)
- Pull the upper and lower headbands over the back of the head whilst positioning the respirator on your face. (Fig 2a)
- Position the upper strap at the crown of your head and the lower strap at the neck. Ensure they are not twisted and are sitting flat. (Fig 3a)
- Pull both ends of the head strap at the same time to evenly tighten the respirator against the nose to a comfortable position (Fig 4a), then do the same for the neck strap to pull against the chin (Fig 5a). Adjust each strap and the position of the respirator on the face accordingly to achieve a sealed face with no air leakage between the respirator and face. (CARE SHOULD BE TAKEN NOT TO OVERTIGHTEN THE STRAPS AS THIS COULD CAUSE THE MASK TO COLLAPSE. WHILST THIS WILL NOT NEGATIVELY IMPACT THE PERFORMANCE, EXCESSIVE STRAP TENSION IS NOT NECESSARY TO ACHIEVE THE BEST FACE FIT).
- Perform a positive pressure face check by exhaling deeply and feeling for air leakage around the nose and face to ensure an air tight seal has been achieved. Perform this each time the respirator is re-applied and repeat steps C – G if a leak is found. (Fig 6a)

USE

NAFAShield 99 FFP3 V R D Safety Respirator is suitable for multipleshift use, however it is recommended that the respirator is inspected before each use. Replace the respirator as frequently as the application requires or immediately if breathing resistance increases during use.

CLEANING AND DISINFECTION

Using a soft cloth dampened with a mild solution of neutral detergent and water, wipe the respirator body and harness assembly avoiding contact with the filtration material.

The filtration material is not washable.

RE-USING THE HARNESS – (MULTI-PACK ONLY)

Only one harness assembly is supplied in a multi-pack, therefore it must be re-used with the additional respirator bodies supplied. To transfer the harness onto a new respirator body, follow the steps described below:

- Looking at the outside of the respirator, pull out the securing plug. (Fig 1b)
- Pull bottom of harness away from the respirator body and slide out of the top nose fixture. (Fig 2b)
- Offer new respirator body to the top of harness and slide under the top nose fixture. (Fig 3b)
- Pull the bottom of the harness down over the exhalation valve ensuring that it sits behind the highlighted securing lugs on both sides. (Fig 4b)
- Push the securing plug back into place the harness assembly onto the respirator body. Ensure the plug is fully inserted. Pull the plug into position from the inside of the mask if necessary. (Fig 5b)

STORAGE

- NAFAShield 99 FFP3 V R D Safety Respirator should be stored in a dry and uncontaminated environment between temperatures of -10 to 40 degrees C.
- NAFAShield 99 FFP3 V R D Safety Respirator should not be stored for extended periods in direct sunlight or beyond the specified shelf life.
- If stored under conditions other than as recommended by the manufacturer the shelf life of the mask may be affected.

ADDITIONAL WARNINGS & LIMITATIONS

DO NOT use the NAFAShield 99 FFP3 V R D Safety Respirator if:

- the air quality, the atmosphere or the type of hazard are unknown.
- the atmosphere is oxygen (<19.5%).
- the atmosphere is oxygen enriched.
- the gas or vapour concentrations are unknown or exceed the value of the class of the mask.
- the atmosphere is explosive.
- the intended user has facial hair protruding beyond the border.
- the mask has exceeded its expiry date, which is displayed on the mask pack label.

ALWAYS:

- Check the and effective seal of the facepiece prior to use – Note: it is unlikely that the requirement for leakage will be achieved if facial hair passes under the face seal.
- Check the respirator for damage and replace if in doubt.
- Ensure that the intended application and atmosphere has been assessed by a competent person and meets the application limitations.
- Check the expiry date, the detail is displayed on the pack label.
- *NOTE* - Do not wash the filters.

FITTING INSTRUCTION DIAGRAM GAMBARAJAH CARA PEMAKAIAN



DIAGRAM CHANGING THE HARNESS GAMBARAJAH PENUKARAN TALI PEMEGANG

PRODUCT DISCLAIMER

NAFAShield 99 FFP3 V R D product offers are correct as at issuance of this document or in any of our related publications. The technical information presented herein is correct and best to our knowledge. However, the manufacturer reserves the right to make changes where applicable without any form of prior notice. The NAFAShield 99 FFP3 V R D product performance as claimed in this document or in any of our publications, is based on testing under controlled conditions according to the testing regimentation set by and strictly complying to the stated test-standards. The manufacturer and its authorized agents are not responsible for variation in the product performance when applying the device under non-standard non-controlled conditions resulting in product performance differ than those stated herein. Users are strongly advised to be alert on the conditions they are exposed to and any application of this device under those conditions are at their own risks. Users must closely follow all recommended cautionary measures by the relevant authorities and specialists.